

THE THERMOCOUPLE

Clean The Bowl



Deep Clean

750 gms Water
250 gms Vinegar
2 tbs Bi-Carb Soda

15 mins | Varoma | Speed 4



SUE & CHRIS ONSLOW

THE THERMOCOUPLE

Clean The Bowl



Shine & Deodorise

3-4 Empty Egg Shells
1 Lemon Cut Up
10 Seconds | Speed 8
Scrape down the side & lid
Add 500 gms Water
30 Mins | 50 C | Speed 2



SUE & CHRIS ONSLOW

THE THERMOCOUPLE

Quick Clean



Easy Quick Clean of the bowl takes literally 3 seconds

1000 gms Water
Quick spray/1 drop dishwashing liquid

Turbo x 3 then rinse



SUE & CHRIS ONSLOW

THE THERMOCOUPLE

Erase Odours



The Bowl can end up smelly as we all know. This is a great way to refresh the smell

1000 gms Water
Lemon Juice
2 tsp Bi-Carb Soda

15 mins | Varoma | Speed 3



SUE & CHRIS ONSLOW

THE THERMOCOUPLE

Restore The Lid



Got a Yellow Lid of your TM Bowl?
First things first
DON'T PANIC

Wash the lid and take it outside
Place it in the direct sun for about 20 mins - that is usually long enough to bring it back to its original colour and set your heart back to normal pace :)



SUE & CHRIS ONSLOW

THE THERMOCOUPLE

Sterilise



Baby Bottles

You can sterilise your babies bottles, dummies, formula travel containers and more

Add 500 ml Water to TM Bowl
Place items in the Varoma Tray
Secure the lid
15 mins | Varoma | Speed 2



SUE & CHRIS ONSLOW